


A4GNT Antibody (R35144)

Catalog No.	Formulation	Size
R35144-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51146
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This A4GNT antibody is available for research use only.



Western blot testing of K562 lysate with A4GNT antibody at 1ug/ml. Predicted molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Alpha-1,4-N-acetylglucosaminyltransferase

Application Notes

Optimal dilution of the A4GNT antibody should be determined by the researcher.

Immunogen

Amino acids PRTYRDLIKPEGS were used as the immunogen for this A4GNT antibody.

Storage

Aliquot and store the A4GNT antibody at -20oC.