


## AAAS Antibody / Adracalin (F53900)

Catalog No.	Formulation	Size
F53900-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53900-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q9NRG9
<b>Applications</b>	Western Blot : 1:1000-2000
<b>Limitations</b>	This AAAS antibody is available for research use only.



Western blot testing of human 1) 293T/17 and 2) HeLa cell lysate with AAAS antibody at 1:2000. Predicted molecular weight: 60 kDa.

## Description

Achalasia adrenocortical insufficiency alacrimia plays a role in the normal development of the peripheral and central nervous system.

## Application Notes

Titration of the AAAS antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 395-425 from the human protein was used as the immunogen for the AAAS antibody.

## Storage

Aliquot the AAAS antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.