


AADAT Antibody (R33618)

Catalog No.	Formulation	Size
R33618-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51166
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AADAT antibody is available for research use only.



IHC testing of FFPE human kidney with AADAT antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Amino adipate aminotransferase

Application Notes

Optimal dilution of the AADAT antibody should be determined by the researcher.

Immunogen

Amino acids DRVIDFYSNQKD were used as the immunogen for this AADAT antibody.

Storage

Aliquot and store the AADAT antibody at -20oC.