


## AAK1 Antibody (R34135)

Catalog No.	Formulation	Size
R34135-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22848
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This AAK1 antibody is available for research use only.



Western blot testing of human amygdala lysate with AAK1 antibody at 1ug/ml. Predicted molecular weight: ~104kDa.

## Description

Additional name(s) for this target protein: AP2 associated kinase 1

## Application Notes

Optimal dilution of the AAK1 antibody should be determined by the researcher.

## Immunogen

Amino acids GSPRTSQQNVYNPSE were used as the immunogen for this AAK1 antibody.

## Storage

Aliquot and store the AAK1 antibody at -20oC.