


## ABCA9 Antibody (R35449)


Catalog No.	Formulation	Size
R35449-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10350
Applications	Western Blot : 0.1-0.3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ABCA9 antibody is available for research use only.



Western blot testing of cerebellum lysate with ABCA9 antibody at 0.1ug/ml. Predicted molecular weight: ~184 kDa.



IHC testing of FFPE human testis tissue with ABCA9 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family A (ABC1), member 9

## Application Notes

Optimal dilution of the ABCA9 antibody should be determined by the researcher.

## Immunogen

Amino acids QRVVQELEMENIQD were used as the immunogen for this ABCA9 antibody.

## Storage

Aliquot and store the ABCA9 antibody at -20oC.