


ABCA9 Antibody (R35449)


Catalog No.	Formulation	Size
R35449-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10350
Applications	Western Blot : 0.1-0.3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ABCA9 antibody is available for research use only.



Western blot testing of cerebellum lysate with ABCA9 antibody at 0.1ug/ml. Predicted molecular weight: ~184 kDa.



IHC testing of FFPE human testis tissue with ABCA9 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family A (ABC1), member 9

Application Notes

Optimal dilution of the ABCA9 antibody should be determined by the researcher.

Immunogen

Amino acids QRVVQELEMENIQD were used as the immunogen for this ABCA9 antibody.

Storage

Aliquot and store the ABCA9 antibody at -20oC.