


ABCB5 Antibody (R35491)

Catalog No.	Formulation	Size
R35491-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	340273
Applications	Western Blot : 0.3-1ug/ml Flow Cytometry : suitable ELISA (peptide) LOD : 1:32000
Limitations	This ABCB5 antibody is available for research use only.



Western blot testing of A431 lysate with ABCB5 antibody at 0.3ug/ml. Predicted molecular weight ~90kDa, observed here at ~80kDa.

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family B (MDR/TAP), member 5

Application Notes

Optimal dilution of the ABCB5 antibody should be determined by the researcher.

Immunogen

Amino acids QTQHRNTSKKAQ were used as the immunogen for this ABCB5 antibody.

Storage

Aliquot and store the ABCB5 antibody at -20oC.