


ABCB5 Antibody (R35491)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35491-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 340273 |
| Applications | Western Blot : 0.3-1ug/ml Flow Cytometry : suitable ELISA (peptide) LOD : 1:32000 |
| Limitations | This ABCB5 antibody is available for research use only. |



Western blot testing of A431 lysate with ABCB5 antibody at 0.3ug/ml. Predicted molecular weight ~90kDa, observed here at ~80kDa.

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family B (MDR/TAP), member 5

Application Notes

Optimal dilution of the ABCB5 antibody should be determined by the researcher.

Immunogen

Amino acids QTQHRNTSKKAQ were used as the immunogen for this ABCB5 antibody.

Storage

Aliquot and store the ABCB5 antibody at -20oC.