


ABCB9 Antibody (R34095)

Catalog No.	Formulation	Size
R34095-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23457
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ABCB9 antibody is available for research use only.



IHC testing of FFPE human testis with ABCB9 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family B (MDR/TAP), member 9

Application Notes

Optimal dilution of the ABCB9 antibody should be determined by the researcher.

Immunogen

Amino acids GHNEPVANGSHKA were used as the immunogen for this ABCB9 antibody.

Storage

Aliquot and store the ABCB9 antibody at -20oC.