


ABCC4 Antibody / MRP4 (R35530)

Catalog No.	Formulation	Size
R35530-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10257
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ABCC4 antibody is available for research use only.



Western blot testing of human prostate lysate with ABCC4 antibody at 0.3ug/ml.
Predicted molecular weight: 150 kDa (NP_005836.2), ~97kDa (NP_001098985.1).

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family C (CFTR/MRP), member 4

Application Notes

Optimal dilution of the ABCC4 antibody should be determined by the researcher.

Immunogen

Amino acids RAENDAQKPSLTR were used as the immunogen for this ABCC4 antibody.

Storage

Aliquot and store the ABCC4 antibody at -20oC.