


## ABCD1 Antibody (R33735)

Catalog No.	Formulation	Size
R33735-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	215
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ABCD1 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human ABCD1 and probed with ABCD1 antibody at 1ug/ml (mock transfection in lane 1). The predicted molecular weight is ~83kDa which is observed in transfected lysate but human brain lysate shows a 50kDa band.

## Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family D (ALD), member 1

## Application Notes

Optimal dilution of the ABCD1 antibody should be determined by the researcher.

## Immunogen

Amino acids EDMQRKGYSEQD were used as the immunogen for this ABCD1 antibody.

## Storage

Aliquot and store the ABCD1 antibody at -20oC.