


ABCD2 Antibody (R35069)

Catalog No.	Formulation	Size
R35069-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	225
Localization	Cytoplasmic
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ABCD2 antibody is available for research use only.



Western blot testing of human liver lysate with ABCD2 antibody at 1ug/ml. Predicted molecular weight: ~83 kDa.

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family D (ALD), member 2

Application Notes

Optimal dilution of the ABCD2 antibody should be determined by the researcher.

Immunogen

Amino acids PLSDTLAIKGKVID were used as the immunogen for this ABCD2 antibody.

Storage

Aliquot and store the ABCD2 antibody at -20oC.