


ABCD3 Antibody Biotin Conjugate (R34973BTN)

Catalog No.	Formulation	Size
R34973BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5825
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:16000
Limitations	This ABCD3 antibody is available for research use only.



Western blot testing of human kidney lysate with biotinylated ABCD3 antibody at 1ug/ml.
Predicted molecular weight: ~75kDa.

Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family D

Application Notes

Optimal dilution of the ABCD3 antibody should be determined by the researcher.

1. This ABCD3 antibody is specific for isoform 1.

Immunogen

Amino acids PDGREDQKRKGISD were used as the immunogen for this ABCD3 antibody.

Storage

Aliquot and store the ABCD3 antibody at -20oC.