


## Abeta Antibody (R33372)

Catalog No.	Formulation	Size
R33372-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11820
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Abeta antibody is available for research use only.



Western blot testing of mouse brain lysate with Abeta antibody at 0.03ug/ml. Predicted molecular weight 79~120kDa depending on glycosylation level.

## Description

Additional name(s) for this target protein: Amyloid beta precursor protein; App

## Application Notes

Optimal dilution of the Abeta antibody should be determined by the researcher.

## Immunogen

Amino acids DAEFGHDSGFVVRHQK were used as the immunogen for this Abeta antibody.

## Storage

Aliquot and store the Abeta antibody at -20oC.