


Abeta Antibody (R33372)

Catalog No.	Formulation	Size
R33372-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11820
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Abeta antibody is available for research use only.



Western blot testing of mouse brain lysate with Abeta antibody at 0.03ug/ml. Predicted molecular weight 79~120kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Amyloid beta precursor protein; App

Application Notes

Optimal dilution of the Abeta antibody should be determined by the researcher.

Immunogen

Amino acids DAEFGHDSGFVVRHQK were used as the immunogen for this Abeta antibody.

Storage

Aliquot and store the Abeta antibody at -20oC.