


ABHD12 Antibody (R35584)

Catalog No.	Formulation	Size
R35584-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26090
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ABHD12 antibody is available for research use only.



Western blot testing of mouse brain lysate with ABHD12 antibody at 0.3ug/ml. Predicted molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Abhydrolase domain containing 12

Application Notes

Optimal dilution of the ABHD12 antibody should be determined by the researcher.

1. This ABHD12 antibody is specific for isoform a.

Immunogen

Amino acids REFLGKSEPEHQH were used as the immunogen for this ABHD12 antibody.

Storage

Aliquot and store the ABHD12 antibody at -20oC.