


ABHD14B Antibody (R34246)

Catalog No.	Formulation	Size
R34246-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84836
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ABHD14B antibody is available for research use only.



Western blot testing of human [A], mouse [B], rat [C] liver lysate with ABHD14B antibody at 0.3ug/ml. Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: Abhydrolase domain containing 14B

Application Notes

Optimal dilution of the ABHD14B antibody should be determined by the researcher.

Immunogen

Amino acids HPCYLDKPEEWHT were used as the immunogen for this ABHD14B antibody. No cross-reaction with ABHD14A is expected.

Storage

Aliquot and store the ABHD14B antibody at -20oC.