


ABHD3 Antibody (F43291)

Catalog No.	Formulation	Size
F43291-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43291-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8WU67
Applications	Western Blot : 1:1000
Limitations	This ABHD3 antibody is available for research use only.



Western blot analysis of lysate from human SH-SY5Y cell line using ABHD3 antibody at 1:1000 for each lane. Predicted molecular weight ~46 kDa.

Description

Abhydrolase domain-containing protein 3 contains an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined.

Application Notes

Titration of the ABHD3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 282-310 from the human protein was used as the immunogen for this ABHD3 antibody.

Storage

Aliquot the ABHD3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.