


## ABHD8 Antibody (F41117)

| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F41117-0.4ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml  |
| F41117-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

[Bulk quote request](#)

|                             |                                                         |
|-----------------------------|---------------------------------------------------------|
| <b>Availability</b>         | 1-3 business days                                       |
| <b>Species Reactivity</b>   | Human                                                   |
| <b>Predicted Reactivity</b> | Primate                                                 |
| <b>Format</b>               | Antigen affinity purified                               |
| <b>Host</b>                 | Rabbit                                                  |
| <b>Clonality</b>            | Polyclonal (rabbit origin)                              |
| <b>Isotype</b>              | Rabbit Ig                                               |
| <b>Purity</b>               | Antigen affinity                                        |
| <b>UniProt</b>              | Q96I13                                                  |
| <b>Applications</b>         | Western Blot : 1:1000                                   |
| <b>Limitations</b>          | This ABHD8 antibody is available for research use only. |



## Description

This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. Abhydrolase domain-containing protein 8 contains alpha/beta hydrolase fold and secretory lipase domains.

## Application Notes

Titration of the ABHD8 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 410-439 from the human protein was used as the immunogen for this ABHD8 antibody.

## Storage

Aliquot the ABHD8 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.