

ABHD8 Antibody (F41117)

Catalog No.	Formulation	Size
F41117-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41117-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q96I13
Applications	Western Blot : 1:1000
Limitations	This ABHD8 antibody is available for research use only.

95
72
55
36
28

ABHD8 antibody western blot analysis in human NCI-H292 lysate. Predicted molecular weight ~47 kDa.

Description

This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. Abhydrolase domain-containing protein 8 contains alpha/beta hydrolase fold and secretory lipase domains.

Application Notes

Titration of the ABHD8 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 410-439 from the human protein was used as the immunogen for this ABHD8 antibody.

Storage

Aliquot the ABHD8 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.