


ABI3BP Antibody (F46997)

Catalog No.	Formulation	Size
F46997-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F46997-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q7Z7G0
Applications	Western Blot : 1:1000
Limitations	This ABI3BP antibody is available for research use only.



Description

ABI gene family member 3-binding protein is involved in the regulation of mesenchymal stem cells.

Application Notes

Titration of the ABI3BP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 648-675 from the human protein was used as the immunogen for this ABI3BP antibody.

Storage

Aliquot the ABI3BP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.