


## ACACB Antibody / Acetyl-coenzyme A carboxylase beta (R35351)

Catalog No.	Formulation	Size
R35351-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	32
<b>Applications</b>	Enzyme ImmunoAssay : used as capture for recombinant protein (coat at 1-2ug/ml) ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This ACACB antibody is available for research use only.



Sandwich ELISA with ACACB antibody at 1.5ug/ml as the detect

### Description

Additional name(s) for this target protein: Acetyl-CoA carboxylase 2

### Application Notes

Optimal dilution of the ACACB antibody should be determined by the researcher.

### Immunogen

Amino acids QLGEPLSDKDRKD were used as the immunogen for this ACACB antibody.

## **Storage**

Aliquot and store the ACACB antibody at -20oC.