


ACADVL Antibody (R35772)

Catalog No.	Formulation	Size
R35772-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	37
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ACADVL antibody is available for research use only.



Western blot testing of human heart lysate with ACADVL antibody at 1ug/ml. Predicted molecular weight: ~68 kDa.

Description

Additional name(s) for this target protein: Acyl-Coenzyme A dehydrogenase, very long chain; LCACD; VLCAD

Application Notes

Optimal dilution of the ACADVL antibody should be determined by the researcher.

Immunogen

Amino acids RSLSEGHPTAQHE were used as the immunogen for this ACADVL antibody.

Storage

Aliquot and store the ACADVL antibody at -20oC.