


ACAP1 Antibody (R33838)

Catalog No.	Formulation	Size
R33838-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9744
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ACAP1 antibody is available for research use only.



Western blot testing of Jurkat lysate with ACAP1 antibody at 0.5ug/ml. Predicted molecular weight: ~81kDa.

Description

Additional name(s) for this target protein: Arf GAP with coiled coil, ANK repeat and PH domains 1; CENTB1; centaurin, beta 1

Application Notes

Optimal dilution of the ACAP1 antibody should be determined by the researcher.

Immunogen

Amino acids EKLSRRSHDLHTL were used as the immunogen for this ACAP1 antibody.

Storage

Aliquot and store the ACAP1 antibody at -20oC.