


ACAT1 Antibody Biotin Conjugate (R33414BTN)

Catalog No.	Formulation	Size
R33414BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	38
Applications	Western Blot: 0.1-0.3ug/ml ELISA (peptide) LOD: 1:2000
Limitations	This ACAT1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Acetyl-Coenzyme A acetyltransferase 1

Application Notes

Optimal dilution of the ACAT1 antibody should be determined by the researcher.

Immunogen

Amino acids DEEYKRVDFSKVPK were used as the immunogen for this ACAT1 antibody. **Storage** Aliquot and store the ACAT1 antibody at -20oC.