


## ACAT1 Antibody Biotin Conjugate (R34636BTN)

Catalog No.	Formulation	Size
R34636BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	38
Applications	Western Blot : 0.1-0.3ug/ml Immunofluorescence : suitable ELISA (peptide) LOD : 1:8000
Limitations	This ACAT1 antibody is available for research use only.



Western blot testing of mouse liver lysate with biotinylated ACAT1 antibody at 0.1ug/ml.  
Predicted molecular weight: ~45kDa.

### Description

Additional name(s) for this target protein: Acetyl-Coenzyme A acetyltransferase 1

### Application Notes

Optimal dilution of the ACAT1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids KRVDFSKVPKLKT were used as the immunogen for this ACAT1 antibody.

## **Storage**

Aliquot and store the ACAT1 antibody at -20oC.