


ACAT1 Antibody (R33414)

Catalog No.	Formulation	Size
R33414-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	38
Applications	Western Blot : 0.1-0.3ug/ml Immunofluorescence : suitable ELISA (peptide) LOD : 1:2000
Limitations	This ACAT1 antibody is available for research use only.



Western blot testing of mouse [A] and rat [B] liver lysate with ACAT1 antibody at 0.1ug/ml. Predicted molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Acetyl-Coenzyme A acetyltransferase 1

Application Notes

Optimal dilution of the ACAT1 antibody should be determined by the researcher.

Immunogen

Amino acids DEEYKRVDFSKVPK were used as the immunogen for this ACAT1 antibody.

Storage

Aliquot and store the ACAT1 antibody at -20oC.