


ACOX2 Antibody (R34261)

Catalog No.	Formulation	Size
R34261-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8309
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ACOX2 antibody is available for research use only.



Western blot testing of human liver lysate with ACOX2 antibody at 0.1ug/ml. Molecular weight: the expected ~75kDa band and the additional ~150kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Acyl-Coenzyme A oxidase 2

Application Notes

Optimal dilution of the ACOX2 antibody should be determined by the researcher.

Immunogen

Amino acids HQSRLRPSDPEAK were used as the immunogen for this ACOX2 antibody.

Storage

Aliquot and store the ACOX2 antibody at -20oC.