


ACSS3 Antibody (F51213)

Catalog No.	Formulation	Size
F51213-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51213-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9H6R3
Applications	Western Blot : 1:500-1:1000
Limitations	This ACSS3 antibody is available for research use only.



Description

Acyl-CoA synthetase short-chain family member 3 activates acetate so that it can be used for lipid synthesis or for energy generation.

Application Notes

Titration of the ACSS3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 224-251 from the human protein was used as the immunogen for this ACSS3 antibody.

Storage

Aliquot the ACSS3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.