


ACT1 Antibody (R34983)

Catalog No.	Formulation	Size
R34983-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O43734
Gene ID	10758
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This ACT1 antibody is available for research use only.



Western blot testing of MCF7 lysate with ACT1 antibody at 2ug/ml. Predicted molecular weight: ~64kDa. (1)

Description

Additional name(s) for this target protein: TRAF3IP2, TRAF3 interacting protein 2, CIKS, ; TRAF3IP2

Application Notes

Optimal dilution of the ACT1 antibody should be determined by the researcher.

1. This Ab will detect isoforms 1 (574 AA) and 2 (565 AA) per UniProt O43734.

Immunogen

Amino acids PDTGHDSKSDQ were used as the immunogen for this ACT1 antibody.

Storage

Aliquot and store the ACT1 antibody at -20°C.