


Actin-like 7B Antibody (R35356)

Catalog No.	Formulation	Size
R35356-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10880
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Actin-like 7B antibody is available for research use only.



Western blot testing of human testis lysate with Actin-like 7B antibody at 0.03ug/ml.
Predicted molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: ACTL7B

Application Notes

Optimal dilution of the Actin-like 7B antibody should be determined by the researcher.

Immunogen

Amino acids QLKMKPRKVHKIK were used as the immunogen for this Actin-like 7B antibody.

Storage

Aliquot and store the Actin-like 7B antibody at -20oC.