


ACTL6A Antibody / ACTL6B (R36356)

Catalog No.	Formulation	Size
R36356-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	86
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ACTL6A antibody is available for research use only.



Western blot testing of HeLa lysate with ACTL6A antibody at 1ug/ml. Expected molecular weight: 47-53 kDa.

Description

Additional name(s) for this target protein: Actin-like 6A; actin-like 6B; 53 kDa BRG1-associated factor A; 53 kDa BRG1-associated factor B

Application Notes

Optimal dilution of the ACTL6A antibody should be determined by the researcher.

Immunogen

Amino acids YEEGGKQCVERKCP were used as the immunogen for this ACTL6A antibody. This Ab will also detect the highly homologous ACTL6B.

Storage

Aliquot and store the ACTL6A antibody at -20oC.