


## Actl7b Antibody (R33419)

Catalog No.	Formulation	Size
R33419-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11471
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Actl7b antibody is available for research use only.



Western blot testing of rat testis lysate with Actl7b antibody at 0.5ug/ml. The expected ~45kDa band and the additional ~30kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Actin-like 7B

## Application Notes

Optimal dilution of the Actl7b antibody should be determined by the researcher.

1. This Actl7b antibody does not react with Actl7a.

## Immunogen

Amino acids DELHVDYELPDGK were used as the immunogen for this Actl7b antibody.

## Storage

Aliquot and store the Actl7b antibody at -20oC.