


## ACTR1A Antibody (R36354)

Catalog No.	Formulation	Size
R36354-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10121
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This ACTR1A antibody is available for research use only.



Western blot testing of human kidney lysate with ACTR1A antibody at 2ug/ml. Predicted molecular weight: ~43kDa.

## Description

Additional name(s) for this target protein: ARP1 actin-related protein 1A

## Application Notes

Optimal dilution of the ACTR1A antibody should be determined by the researcher.

## Immunogen

Amino acids YEEDGARSIHRKT were used as the immunogen for this ACTR1A antibody.

## Storage

Aliquot and store the ACTR1A antibody at -20oC.