


## ADH5 Antibody (R34494)

Catalog No.	Formulation	Size
R34494-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	128
Applications	Western Blot : 0.5-1.5ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:4000
Limitations	This ADH5 antibody is available for research use only.



Western blot testing of human testis lysate with ADH5 antibody at 0.5ug/ml. Predicted molecular weight: ~40kDa.

## Description

Additional name(s) for this target protein: Alcohol dehydrogenase 5

## Application Notes

Optimal dilution of the ADH5 antibody should be determined by the researcher.

## **Immunogen**

Amino acids KKIKVDEFVTHN were used as the immunogen for this ADH5 antibody.

## **Storage**

Aliquot and store the ADH5 antibody at -20oC.