


## ADIPOQ Antibody / Adiponectin [clone 28A36] (RQ8372)

Catalog No.	Formulation	Size
RQ8372	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28A36
Purity	Affinity chromatography
UniProt	Q15848
Applications	Western Blot : 1:500-1:2000
Limitations	This ADIPOQ antibody is available for research use only.



Western blot testing of human plasma cell lysate with ADIPOQ antibody. Predicted molecular weight ~30 kDa.

### Description

Adiponectin is an important adipokine involved in the control of fat metabolism and insulin sensitivity, with direct anti-diabetic, anti-atherogenic and anti-inflammatory activities. [UniProt]

### Application Notes

Optimal dilution of the ADIPOQ antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide derived from human Adiponectin was used as the immunogen for the ADIPOQ antibody.

## **Storage**

Store the ADIPOQ antibody at -20oC.