


ADRA1A Antibody (R35154)

Catalog No.	Formulation	Size
R35154-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	148
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ADRA1A antibody is available for research use only.



Western blot testing of human prostate lysate with ADRA1A antibody at 1ug/ml.
Predicted molecular weight: ~52kDa.

Description

Additional name(s) for this target protein: Alpha adrenergic receptor 1a

Application Notes

Optimal dilution of the ADRA1A antibody should be determined by the researcher.

1. This ADRA1A antibody is specific for isoform 1.

Immunogen

Amino acids PSLDKNHQVPTIK were used as the immunogen for this ADRA1A antibody.

Storage

Aliquot and store the ADRA1A antibody at -20oC.