


ADRA1B Antibody (R36077)


Catalog No.	Formulation	Size
R36077-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	147
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ADRA1B antibody is available for research use only.



Western blot testing of A) human, B) mouse and C) rat heart lysate with ADRA1B antibody at 0.1ug/ml. Predicted molecular weight: ~57kDa.



Western blot of pig heart lysate with ADRA1B antibody at 0.1ug/ml. Predicted molecular weight: ~57kDa.

Description

Additional name(s) for this target protein: Alpha-1B adrenoceptor

Application Notes

Optimal dilution of the ADRA1B antibody should be determined by the researcher.

Immunogen

Amino acids SSTKAKGHNPRSS were used as the immunogen for this ADRA1B antibody.

Storage

Aliquot and store the ADRA1B antibody at -20oC.