


## ADRA1B Antibody (R36077)


Catalog No.	Formulation	Size
R36077-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat, Pig
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	147
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This ADRA1B antibody is available for research use only.



Western blot testing of A) human, B) mouse and C) rat heart lysate with ADRA1B antibody at 0.1ug/ml. Predicted molecular weight: ~57kDa.



Western blot of pig heart lysate with ADRA1B antibody at 0.1ug/ml. Predicted molecular weight: ~57kDa.

## Description

Additional name(s) for this target protein: Alpha-1B adrenoceptor

## Application Notes

Optimal dilution of the ADRA1B antibody should be determined by the researcher.

## Immunogen

Amino acids SSTKAKGHNPRSS were used as the immunogen for this ADRA1B antibody.

## Storage

Aliquot and store the ADRA1B antibody at -20oC.