


## ADRB1 Antibody (R33987)

Catalog No.	Formulation	Size
R33987-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	153
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ADRB1 antibody is available for research use only.



Western blot testing of human heart lysate with ADRB1 antibody at 0.5ug/ml. Predicted molecular weight: ~51 kDa.

## Description

Additional name(s) for this target protein: Beta-1-adrenergic receptor

## Application Notes

Optimal dilution of the ADRB1 antibody should be determined by the researcher.

## Immunogen

Amino acids ESDEARRCYNDPK were used as the immunogen for this ADRB1 antibody.

## **Storage**

Aliquot and store the ADRB1 antibody at -20oC.