

ADRB1 Antibody (R33987)

Catalog No.	Formulation	Size
R33987-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	153
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ADRB1 antibody is available for research use only.



Western blot testing of human heart lysate with ADRB1 antibody at 0.5ug/ml. Predicted molecular weight: ~51 kDa.

Description

Additional name(s) for this target protein: Beta-1-adrenergic receptor

Application Notes

Optimal dilution of the ADRB1 antibody should be determined by the researcher.

Immunogen

Amino acids ESDEARRCYNDPK were used as the immunogen for this ADRB1 antibody.

Storage

Aliquot and store the ADRB1 antibody at -20oC.