


## ADRB3 Antibody (R35019)

Catalog No.	Formulation	Size
R35019-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	155
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ADRB3 antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with ADRB3 antibody at 0.03ug/ml.  
Predicted molecular weight: ~44kDa.

## Description

Additional name(s) for this target protein: Adrenergic beta-3 receptor antibody, BETA3AR, Beta-3 Adrenergic Receptor

## Application Notes

Optimal dilution of the ADRB3 antibody should be determined by the researcher.

## Immunogen

Amino acids PHENSSLAPWPD were used as the immunogen for this ADRB3 antibody.

## Storage

Aliquot and store the ADRB3 antibody at -20oC.