


ADRB3 Antibody (R35019)

Catalog No.	Formulation	Size
R35019-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	155
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ADRB3 antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with ADRB3 antibody at 0.03ug/ml.
Predicted molecular weight: ~44kDa.

Description

Additional name(s) for this target protein: Adrenergic beta-3 receptor antibody, BETA3AR, Beta-3 Adrenergic Receptor

Application Notes

Optimal dilution of the ADRB3 antibody should be determined by the researcher.

Immunogen

Amino acids PHENSSLAPWPD were used as the immunogen for this ADRB3 antibody.

Storage

Aliquot and store the ADRB3 antibody at -20oC.