


AGAP3 Antibody (R34881)

Catalog No.	Formulation	Size
R34881-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	116988
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AGAP3 antibody is available for research use only.



IHC testing of FFPE human kidney with AGAP3 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: CENTG3, MRIP1

Application Notes

Optimal dilution of the AGAP3 antibody should be determined by the researcher.

Immunogen

Amino acids NPSAELHRSPSLL were used as the immunogen for this AGAP3 antibody.

Storage

Aliquot and store the AGAP3 antibody at -20oC.