


AGAP3 Antibody (R34881)

Catalog No.	Formulation	Size
R34881-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	116988
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AGAP3 antibody is available for research use only.



IHC testing of FFPE human kidney with AGAP3 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: CENTG3, MRIP1

Application Notes

Optimal dilution of the AGAP3 antibody should be determined by the researcher.

Immunogen

Amino acids NPSAELHRSPSLL were used as the immunogen for this AGAP3 antibody. **Storage** Aliquot and store the AGAP3 antibody at -20oC.