


Agbl5 Antibody [isoform 3] (R33232)

Catalog No.	Formulation	Size
R33232-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	231093
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Agbl5 antibody is available for research use only.



Western blot testing of mouse skin lysate with Agbl5 antibody at 1ug/ml. Predicted molecular weight: ~93kDa.

Description

Additional name(s) for this target protein: ATP/GTP binding protein-like 5

Application Notes

Optimal dilution of the Agbl5 antibody should be determined by the researcher.

1. This Agbl5 antibody is specific to isoform 3 only.

Immunogen

Amino acids CFSKPEEAGSHVE were used as the immunogen for this Agbl5 antibody. (1)

Storage

Aliquot and store the Agbl5 antibody at -20oC.