


AGR2 Antibody [clone GCI-1] (RQ4822)

Catalog No.	Formulation	Size
RQ4822	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	GCI-1
Purity	Affinity purified
UniProt	O95994
Applications	Western Blot : 1:1000-1:2000
Limitations	This AGR2 antibody is available for research use only.



Western blot testing of human MCF7 cell lysate with AGR2 antibody. Expected molecular weight: 17-20 kDa.

Description

Anterior gradient protein 2 plays a role in cell migration, cellular transformation and metastasis and is a p53 inhibitor.

Application Notes

Optimal dilution of the AGR2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human AGR2 was used as the immunogen for the AGR2 antibody.

Storage

Store the AGR2 antibody at -20oC.