


## AGTPBP1 Antibody (R36294)

Catalog No.	Formulation	Size
R36294-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23287
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This AGTPBP1 antibody is available for research use only.



Western blot testing of HeLa lysate with AGTPBP1 antibody at 1ug/ml. Predicted molecular weight: ~138kDa.

### Description

Additional name(s) for this target protein: ATP/GTP binding protein 1; NNA1

### Application Notes

Optimal dilution of the AGTPBP1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids TSPLEYNLPS were used as the immunogen for this AGTPBP1 antibody.

## **Storage**

Aliquot and store the AGTPBP1 antibody at -20oC.