


AGTR1 Antibody (R33808)

Catalog No.	Formulation	Size
R33808-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	185
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AGTR1 antibody is available for research use only.



Western blot testing of mouse liver lysate with AGTR1 antibody at 0.3ug/ml. Observed molecular weight: 41-50 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: AT1, AT1R, Angiotensin II receptor type 1

Application Notes

Optimal dilution of the AGTR1 antibody should be determined by the researcher.

Immunogen

Amino acids EIQKNKPRNDDIFK were used as the immunogen for this AGTR1 antibody.

Storage

Aliquot and store the AGTR1 antibody at -20oC.