


## AGTR1 Antibody (R33808)

Catalog No.	Formulation	Size
R33808-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	185
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AGTR1 antibody is available for research use only.



Western blot testing of mouse liver lysate with AGTR1 antibody at 0.3ug/ml. Observed molecular weight: 41-50 kDa depending on glycosylation level.

## Description

Additional name(s) for this target protein: AT1, AT1R, Angiotensin II receptor type 1

## Application Notes

Optimal dilution of the AGTR1 antibody should be determined by the researcher.

## Immunogen

Amino acids EIQKNKPRNDDIFK were used as the immunogen for this AGTR1 antibody.

## Storage

Aliquot and store the AGTR1 antibody at -20oC.